



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

# **DISTAL FEMUR OSTOTOMY LOCKED PLATE**





**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

# DISTAL FEMUR OSTOTOMY LOCKED PLATE

## DESCRIPTION

<b>Cat No. SS :- S23.DFOLPSS</b>	<b>Cat No. TIT :- S23.DFOLPTT</b>
<b>Suitable For</b>	<b>Fixation of Closed - Wedge Osteotomies Close To The Knee</b>
<b>Screw</b>	<b>5.0 MM Cortex Locked Screw, 4.5 MM Cortical Screw</b>
<b>Hole</b>	<b>4 Hole To 8 Hole</b>



<b>HOLE</b>	<b>CAT. No. S.S :- S23.DFOLPSS</b>	<b>CAT. No. TIT :- S23.DFOLPTT</b>
<b>04 HOLE</b>	<b>S23.SS59280</b>	<b>S23.TIT59280</b>
<b>05 HOLE</b>	<b>S23.SS59281</b>	<b>S23.TIT59281</b>
<b>06 HOLE</b>	<b>S23.SS59282</b>	<b>S23.TIT59282</b>
<b>07 HOLE</b>	<b>S23.SS59283</b>	<b>S23.TIT59283</b>
<b>08 HOLE</b>	<b>S23.SS59284</b>	<b>S23.TIT59284</b>