



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

DORSAL DISTAL RADIUS L- LOCKED PLATE (LEFT & RIGHT)



DORSAL DISTAL RADIUS L- LOCKED PLATE (LEFT & RIGHT)

DESCRIPTION

Cat No. SS :- S23.DDRLLPLSS / S23.DDRLLPRSS	Cat No. TIT :- S23.DDRLLPLTT / S23.DDRLLPRTT
Suitable For	Intra And Extra Articular Fractures And Osteotomies Of Distal Radius
Screw Hole	2.4 MM Cortex Locked Screw, 2.7 MM Cortical Screw 2 x 4 Hole & 2 x 18 Hole



HOLE	LEFT		RIGHT	
	CAT. No. S.S :- S23.DDRLLPLSS	CAT. No. TIT :- S23.DDRLLPLTT	CAT. No. S.S :- S23.DDRLLPRSS	CAT. No. TIT :- S23.DDRLLPRTT
2 X 04 HOLE	S23.SS00548	S23.TIT00548	S23.SS00549	S23.TIT00549
2 X 05 HOLE	S23.SS17377	S23.TIT17377	S23.SS17396	S23.TIT17396
2 X 06 HOLE	S23.SS17378	S23.TIT17378	S23.SS17397	S23.TIT17397
2 X 07 HOLE	S23.SS17379	S23.TIT17379	S23.SS17398	S23.TIT17398
2 X 08 HOLE	S23.SS17380	S23.TIT17380	S23.SS17399	S23.TIT17399
2 X 09 HOLE	S23.SS17381	S23.TIT17381	S23.SS17400	S23.TIT17400
2 X 10 HOLE	S23.SS17382	S23.TIT17382	S23.SS17401	S23.TIT17401
2 X 11 HOLE	S23.SS17383	S23.TIT17383	S23.SS17402	S23.TIT17402
2 X 12 HOLE	S23.SS17384	S23.TIT17384	S23.SS17403	S23.TIT17403
2 X 13 HOLE	S23.SS17385	S23.TIT17385	S23.SS17404	S23.TIT17404
2 X 14 HOLE	S23.SS17386	S23.TIT17386	S23.SS17405	S23.TIT17405
2 X 15 HOLE	S23.SS17387	S23.TIT17387	S23.SS17406	S23.TIT17406
2 X 16 HOLE	S23.SS17388	S23.TIT17388	S23.SS17407	S23.TIT17407
2 X 17 HOLE	S23.SS17389	S23.TIT17389	S23.SS17408	S23.TIT17408
2 X 18 HOLE	S23.SS17390	S23.TIT17390	S23.SS17409	S23.TIT17409