



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

# OLECRANON LOCKED PLATE



## DESCRIPTION

<b>Cat No. SS :- S23.OLPSS</b>	<b>Cat No. TIT :- S23.OLPTT</b>
<b>Suitable For</b>	<b>Complex And Simple Olecranon Fracture</b>
<b>Screw</b>	<b>4.0 MM Cancellous Screw, 4.0 MM Cortex Locked Screw, 3.5 MM Cortical Screw</b>
<b>Hole</b>	<b>3 Hole To 12 Hole</b>



<b>HOLE</b>	<b>CAT. No. S.S :- S23.OLPSS</b>	<b>CAT. No. TIT :- S23.OLPTT</b>
<b>03 HOLE</b>	<b>S23.SS00455</b>	<b>S23.TIT00455</b>
<b>04 HOLE</b>	<b>S23.SS00456</b>	<b>S23.TIT00456</b>
<b>05 HOLE</b>	<b>S23.SS00457</b>	<b>S23.TIT00457</b>
<b>06 HOLE</b>	<b>S23.SS00458</b>	<b>S23.TIT00458</b>
<b>07 HOLE</b>	<b>S23.SS00459</b>	<b>S23.TIT00459</b>
<b>08 HOLE</b>	<b>S23.SS00460</b>	<b>S23.TIT00460</b>
<b>09 HOLE</b>	<b>S23.SS00461</b>	<b>S23.TIT00461</b>
<b>10 HOLE</b>	<b>S23.SS00462</b>	<b>S23.TIT00462</b>
<b>11 HOLE</b>	<b>S23.SS00463</b>	<b>S23.TIT00463</b>
<b>12 HOLE</b>	<b>S23.SS00464</b>	<b>S23.TIT00464</b>