



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**



**Volar Extra Articular Standard  
Locked Plate ( Left & Right )**



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

## DESCRIPTION

<b>Cat No. SS:-S23.VEASLPLSS</b>	<b>Cat No. TIT :- S23.VEASLPRTT</b>
<b>Suitable For</b>	<b>Intra And Extra Articular Fractures And Osteotomies Of Distal Radius</b>
<b>Screw</b>	<b>2.4 MM Cortex Locked Screw, 2.7 MM Cortical Screw</b>
<b>Hole</b>	<b>5 x 4 Hole</b>

	LEFT		RIGHT	
HOLE	CAT. No. S.S :- S23.DDRRALLPLSS	CAT. No. TIT :- S23.DDRRALLPLTT	CAT. No. S.S :- S23.DDRRALLPRSS	CAT. No. TIT :- S23.DDRRALLPRTT
<b>5 X 4 HOLE</b>	<b>S23.SS00551</b>	<b>S23.TIT00551</b>	<b>S23.SS00552</b>	<b>S23.TIT00552</b>

**Volar Extra Articular Standard Locked Plate ( Left & Right )**