



DESCRIPTION

Suitable For	Trochanteric Fracture of Femur
Screw Ø	12.5 MM & 4.5 MM
Profile	19.0 MM X 6.0 MM
Distance Between Hole	16.0 MM
Material	Stainless Steel 316L / 316LVM / TIT.
Hole	2 Hole To 18 Hole (Diff. 1 Hole)
Degree	125,130,135,140
Proximal Screw	LAG SCREW 12.5 MM
Screw Length	25.0 MM To 120.0 MM (Diff.5.0 MM)
Distal Screw	CORTICAL SCREW 4.5 MM
Screw Length	10.0 MM To 54.0 MM (Diff. 2.0 MM)

**Dynamic Compression Plate
L/B Or S/B (R/H,DC/H)**



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

LENGTH	CAT. No. S.S.	CAT. No. TIT.
02 HOLE	S23.SS01297	S23.TIT01297
03 HOLE	S23.SS01298	S23.TIT01298
04 HOLE	S23.SS01299	S23.TIT01299
05 HOLE	S23.SS01300	S23.TIT01300
06 HOLE	S23.SS01301	S23.TIT01301
07 HOLE	S23.SS01302	S23.TIT01302
08 HOLE	S23.SS01303	S23.TIT01303
09 HOLE	S23.SS01304	S23.TIT01304
10 HOLE	S23.SS01305	S23.TIT01305
11 HOLE	S23.SS01306	S23.TIT01306
12 HOLE	S23.SS01307	S23.TIT01307

LENGTH	CAT. No. S.S.	CAT. No. TIT.
13 HOLE	S23.SS01308	S23.TIT01308
14 HOLE	S23.SS01309	S23.TIT01309
15 HOLE	S23.SS01310	S23.TIT01310
16 HOLE	S23.SS01311	S23.TIT01311
17 HOLE	S23.SS01312	S23.TIT01312
18 HOLE	S23.SS01313	S23.TIT01313

Dynamic Compression Plate L/B Or S/B (R/H,DC/H)