



DESCRIPTION

| | |
|------------------------------|---|
| Suitable For | Supracondylar Fracture Distal Femur Fracture |
| Screw Ø | 4.5 MM & 6.5 MM, 12.5 MM |
| Profile | 6.0 MM & 16.0 MM |
| Distance Between Hole | 16.0 MM |
| Material | Stainless Steel 316L / 316LVM / TIT. |
| Hole | 2 Hole To 16 Hole (Diff. 1 Hole) |
| Proximal Screw | Cortical Screw 4.5 MM |
| Screw Length | 10.0 MM To 54.0 MM (Diff. 2.0 MM) |
| Distal Screw | Lag Screw 12.5 MM |
| Screw Length | 25.0 MM To 120.0 MM (DIFF. 5.0 MM) |
| Distal Screw | Cancellous Screw 6.5 MM |
| Screw Length | 25.0 MM To 110.0 MM (Diff. 5.0 MM) |

Dynamic Condylar Plate 95°



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

| LENGTH | CAT. No. S.S. | CAT. No. TIT. |
|----------------|----------------------|----------------------|
| 02 HOLE | S23.SS01227 | S23.TIT01227 |
| 03 HOLE | S23.SS01228 | S23.TIT01228 |
| 04 HOLE | S23.SS01229 | S23.TIT01229 |
| 05 HOLE | S23.SS01230 | S23.TIT01230 |
| 06 HOLE | S23.SS01231 | S23.TIT01231 |
| 07 HOLE | S23.SS01232 | S23.TIT01232 |
| 08 HOLE | S23.SS01233 | S23.TIT01233 |
| 09 HOLE | S23.SS01234 | S23.TIT01234 |
| 10 HOLE | S23.SS01235 | S23.TIT01235 |
| 11 HOLE | S23.SS01236 | S23.TIT01236 |
| 12 HOLE | S23.SS01237 | S23.TIT01237 |
| 13 HOLE | S23.SS01238 | S23.TIT01238 |
| 14 HOLE | S23.SS01239 | S23.TIT01239 |
| 15 HOLE | S23.SS01240 | S23.TIT01240 |
| 16 HOLE | S23.SS01241 | S23.TIT01241 |

Dynamic Condylar Plate 95°