



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

## DESCRIPTION

<b>Thread Diameter</b>	<b>6.5 MM</b>
<b>Core Diameter</b>	<b>3.8 MM</b>
<b>Use Drills Diameter</b>	<b>4.0 MM</b>
<b>For Passing Hole</b>	<b>6.5 MM</b>
<b>Drive</b>	<b>Hex</b>
<b>Material</b>	<b>Stainless Steel 316L / 316LVM / TIT.</b>
<b>Size</b>	<b>6.5 MM</b>
<b>Length</b>	<b>20.0 MM To 100.0 MM (Diff. 5.0 MM)</b>



**Cancellous Locked Screw  
6.5 mm**



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

LENGTH	16MM THREAD		32MM THREAD		FULL THREAD	
	CAT. No. S.S.	CAT. No. TIT.	CAT. No. S.S.	CAT. No. TIT.	CAT. No. S.S.	CAT. No. TIT.
20.0 MM	S23.SS00755	S23.TIT00755	S23.SS00772	S23.TIT00772	S23.SS00789	S23.TIT00789
25.0 MM	S23.SS00756	S23.TIT00756	S23.SS00773	S23.TIT00773	S23.SS00790	S23.TIT00790
30.0 MM	S23.SS00757	S23.TIT00757	S23.SS00774	S23.TIT00774	S23.SS00791	S23.TIT00791
35.0 MM	S23.SS00758	S23.TIT00758	S23.SS00775	S23.TIT00775	S23.SS00792	S23.TIT00792
40.0 MM	S23.SS00759	S23.TIT00759	S23.SS00776	S23.TIT00776	S23.SS00793	S23.TIT00793
45.0 MM	S23.SS00760	S23.TIT00760	S23.SS00777	S23.TIT00777	S23.SS00794	S23.TIT00794
50.0 MM	S23.SS00761	S23.TIT00761	S23.SS00778	S23.TIT00778	S23.SS00795	S23.TIT00795
55.0 MM	S23.SS00762	S23.TIT00762	S23.SS00779	S23.TIT00779	S23.SS00796	S23.TIT00796

**Cancellous Locked Screw 6.5 mm**



**Dksortho & Sharma  
Pharmaceutical  
Pvt. Ltd.**

	16MM THREAD		32MM THREAD		FULL THREAD	
LENGTH	CAT. No. S.S.	CAT. No. TIT.	CAT. No. S.S.	CAT. No. TIT.	CAT. No. S.S.	CAT. No. TIT.
60.0 MM	S23.SS00763	S23.TIT00763	S23.SS00780	S23.TIT00780	S23.SS00797	S23.TIT00797
65.0 MM	S23.SS00764	S23.TIT00764	S23.SS00781	S23.TIT00781	S23.SS00798	S23.TIT00798
70.0 MM	S23.SS00765	S23.TIT00765	S23.SS00782	S23.TIT00782	S23.SS00799	S23.TIT00799
75.0 MM	S23.SS00766	S23.TIT00766	S23.SS00783	S23.TIT00783	S23.SS00800	S23.TIT00800
80.0 MM	S23.SS00767	S23.TIT00767	S23.SS00784	S23.TIT00784	S23.SS00801	S23.TIT00801
85.0 MM	S23.SS00768	S23.TIT00768	S23.SS00785	S23.TIT00785	S23.SS00802	S23.TIT00802
90.0 MM	S23.SS00769	S23.TIT00769	S23.SS00786	S23.TIT00786	S23.SS00803	S23.TIT00803
95.0 MM	S23.SS00770	S23.TIT00770	S23.SS00787	S23.TIT00787	S23.SS00804	S23.TIT00804
100.0 MM	S23.SS00771	S23.TIT00771	S23.SS00788	S23.TIT00788	S23.SS00805	S23.TIT00805

**Cancellous Locked Screw 6.5 mm**