



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

DESCRIPTION

Thread Diameter	4.0 MM
Core Diameter	2.7 MM
Use Drills Diameter	3.0 MM
For Passing Hole	4.0 MM
Drive	Hex Drive
Material	Stainless Steel 316L / 316LVM / Titanium
Size	4.0 MM
Length	20.0 MM To 64.0 MM (Diff. 2.0 MM)



**Cortex Locked Screw
5.0mm (Solid) (Hex Drive)**



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

LENGTH	CAT. No. S.S.	CAT. No. TIT.
20.0 MM	S23.SS00902	S23.TIT00902
22.0 MM	S23.SS00903	S23.TIT00903
24.0 MM	S23.SS00904	S23.TIT00904
26.0 MM	S23.SS00905	S23.TIT00905
28.0 MM	S23.SS00906	S23.TIT00906
30.0 MM	S23.SS00907	S23.TIT00907
32.0 MM	S23.SS00908	S23.TIT00908
34.0 MM	S23.SS00909	S23.TIT00909
36.0 MM	S23.SS00910	S23.TIT00910
38.0 MM	S23.SS00911	S23.TIT00911
40.0 MM	S23.SS00912	S23.TIT00912
42.0 MM	S23.SS00913	S23.TIT00913

LENGTH	CAT. No. S.S.	CAT. No. TIT.
44.0 MM	S23.SS00914	S23.TIT00914
46.0 MM	S23.SS00915	S23.TIT00915
48.0 MM	S23.SS00916	S23.TIT00916
50.0 MM	S23.SS00917	S23.TIT00917
52.0 MM	S23.SS00918	S23.TIT00918
54.0 MM	S23.SS00919	S23.TIT00919
56.0 MM	S23.SS00920	S23.TIT00920
58.0 MM	S23.SS00921	S23.TIT00921
60.0 MM	S23.SS00922	S23.TIT00922
62.0 MM	S23.SS00923	S23.TIT00923
64.0 MM	S23.SS00924	S23.TIT00924

Cortex Locked Screw 5.0mm (Solid) (Hex Drive)