



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**



3D Locked Plate



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

DESCRIPTION

Cat No. SS :- S23.3DLPSS	Cat No. TIT :- S23.3DLPTT
Suitable For	1st MP J Arthrodesis With Or Without Interposition Bone Comfort Lapidus Fusion / Metatarsal Fractures
Screw	2.4 MM Or 3.0 MM Cortex Locked Screw, 2.7 MM Or 3.5 MM Cortical Screw
Hole	Standard

HOLE	CAT. No. S.S :- S23.3DLPSS	CAT. No. TIT :- S23.3DLPTT
LEFT	S23.SS11444	S23.TIT11444
RIGHT	S23.SS11445	S23.TIT11445

3D Locked Plate