



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**



**Elastic Nail
(Flexible Nail)**



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

DESCRIPTION

Cat No. SS :- S23.FNSS	Cat No. SS :- S23.FNTT
Suitable For	Humerus, Tibia, Femur Fracture Use Pediatric Patient
Nail Diameter	1.5 MM, 2.0 MM, 2.5 MM, 3.0 MM, 3.5 MM, 4.0 MM, 4.5 MM, 5.0 MM
Nail Length	44.0 CM Standard

LENGTH	44.0 CM STANDARD	
DIAMETER	CAT. No. S.S.	CAT. No. TIT.
1.5 MM	S23.SS02803	S23.TIT02803
2.0 MM	S23.SS02804	S23.TIT02804
2.5 MM	S23.SS02805	S23.TIT02805
3.0 MM	S23.SS02806	S23.TIT02806
3.5 MM	S23.SS02807	S23.TIT02807
4.0 MM	S23.SS02808	S23.TIT02808
4.5 MM	S23.SS02809	S23.TIT02809
5.0 MM	S23.SS02810	S23.TIT02810