




**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

L - BUTTRESS PLATE (LEFT & RIGHT)



DESCRIPTION

Cat No. SS :- S23.LBPLSS / S23.LBPRSS	Cat No. TIT :- S23.LBPLTT / S23.LBPRTT
Suitable For	Proximal Lateral Tibia Fracture, Lateral Condylar Fracture
Screw	6.5 MM Cancellous Screw, 4.5 MM Cortical Screw
Hole	3 Hole To 12 Hole



HOLE	LEFT		RIGHT	
	CAT. No. S.S :- S23.LBPLSS	CAT. No. TIT :- S23.LBPLTT	CAT. No. S.S :- S23.LBPRSS	CAT. No. TIT :- S23.LBPRTT
03 HOLE	S23.SS01550	S23.TIT01550	S23.SS01560	S23.TIT01560
04 HOLE	S23.SS01551	S23.TIT01551	S23.SS01561	S23.TIT01561
05 HOLE	S23.SS01552	S23.TIT01552	S23.SS01562	S23.TIT01562
06 HOLE	S23.SS01553	S23.TIT01553	S23.SS01563	S23.TIT01563
07 HOLE	S23.SS01554	S23.TIT01554	S23.SS01564	S23.TIT01564
08 HOLE	S23.SS01555	S23.TIT01555	S23.SS01565	S23.TIT01565
09 HOLE	S23.SS01556	S23.TIT01556	S23.SS01566	S23.TIT01566
10 HOLE	S23.SS01557	S23.TIT01557	S23.SS01567	S23.TIT01567
11 HOLE	S23.SS01558	S23.TIT01558	S23.SS01568	S23.TIT01568
12 HOLE	S23.SS01559	S23.TIT01559	S23.SS01569	S23.TIT01569