



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**



**ACETABULAR CUP
(UN - CEMENTED)**



**Dksortho & Sharma
Pharmaceutical
Pvt. Ltd.**

DESCRIPTION

Material	S.S / TIT.
Head Coverage	165°
Coating	S.S - HA
Coating	TIT - POROUS

LENGTH	CAT. No. S.S.	CAT. No. TIT.
40.0 MM	S23.THR4900	S23.THR4900-T
42.0 MM	S23.THR4896	S23.THR4896-T
44.0 MM	S23.THR0391	S23.THR0391-T
46.0 MM	S23.THR0393	S23.THR0393-T
48.0 MM	S23.THR0395	S23.THR0395-T
50.0 MM	S23.THR0397	S23.THR0397-T
52.0 MM	S23.THR0399	S23.THR0399-T
54.0 MM	S23.THR0401	S23.THR0401-T
56.0 MM	S23.THR0403	S23.THR0403-T
58.0 MM	S23.THR0405	S23.THR0405-T
60.0 MM	S23.THR0407	S23.THR0407-T
62.0 MM	S23.THR0409	S23.THR0409-T
64.0 MM	S23.THR0411	S23.THR0411-T
66.0 MM	S23.THR0413	S23.THR0413-T

ACETABULAR CUP (UN - CEMENTED)